



5024 Rose St * Schiller Park IL 60176-1023
800-442-8263 * Fax 800-442-0487
www.gelatininnovations.com
sales@gelatininnovations.com

PECTIN SPECIFICATION

*Viscosity of 1% solution 40-60 cps
Moisture – Less than 10%
Acid Insoluble Ash Less than 1%
Methoxl Content – Min 6.7%
Galacturonic Acid Content – Min74%

MICRO

Total Plate Count – Less than 1000/g
Yeast and Mold – Less than 100/g
Coliforms – Less than 0.3/g
Salmonella – Negative in 25g

METALS

Arsenic Less than 3 ppm
Lead Less than 5 ppm

VISCOSITY TEST METHOD

1% Pectin/Water Solution; Brookfield LVT; 60 rpm; Spindle #1; 25°C.